Contents

I. Product explanation.....	02
General Information	02
Intended use.....	02
Characteristics.....	02
Components.....	02
Specifications (Technical Specifications).....	04
Performance	04
Symbols	04
II. Instructions for use	07
Frequently Used Function	07
Workflow Diagram.....	07
Preparation Before Use	07
Charging Battery (Charging adapter sold separately).....	08
Operation Method.....	08
Measuring Process.....	09
Checking Measured Value (iST scale).....	10
Product Performance Test.	10
Precautions During Use.....	11
Cleaning After Use.....	11
Sterilization After cleaning	11
Storage after sterilization.....	12
Assembling/Disassembling Method	12
Charging battery.....	13
Storage and Maintenance After Use (Instructions to stop operation).....	13
Disposal of Products and Components.....	13
III. Warnings and Safety Precautions	14
Precaution	14
Prohibition.....	14
Warnings.....	14
IV. Statements and tables for EMC	15
V. Statements for FCC	18
VI. Product Warranty Policy	19

I. Product explanation

0 General Information

Manufacturer	DMS Co., Ltd.
Manufacturer address, Contact Information	4dong-#107, 201, 202, 42-10, Taejunggongdan- gil, Wonju-si, Gangwon-do, Republic of Korea A/S Representative Tel. # : (+82) 33-742-3210
Product Name	AnyCheck
Model Name	IMT-100
Rating	3.7 Vdc, 640 mAh (Charging: 5 Vdc, 1-2 A)
Expected Service life	5 years from the manufactured date
Free Warranty Period	1 year from date of purchase (Even during the free warranty period, costs can be incurred for damages caused by consumers and failures caused by insufficient management)
Software Name	IMT-SW
Software Version	IMT-1.1.X ¹⁾

1) X : This notation is a two-digit number to manage the minor changes of software due to internal bugs or stabilization. The number is marked such as "X" just for our internal management because it is frequently changed. That is not related to the performance and function changes of the product and there is no permission problem.


1 Intended use

IMT-100 is a stability measuring device that uses tapping method, which measures the time of the Tapping Rod of the device contacting the implant for dentistry and dental implantology or oral cavity, or abutment, to measure the stability of the placed implant. The result of the measurement is displayed in iST scale (Implant Stability Test), which is the numerical representation of the Implant Stability Test (iST) value, and higher scale represents greater implant stability or lower implant mobility. And the intended user is dentist and dental nurse.

2 Characteristics

- ① Convenient usage due to small and light design
- ② Easy to use due to simple controlling method
- ③ Fast measuring time less than 3 seconds
- ④ Safe measurement due to slight tapping technique and auto safety system (tapping stops after tapping twice when Osseointegration is unstable, otherwise taps up to 6 times.)
- ⑤ Available to check possibility of implant failure at in earlier stages.

3 Components

I. Product explanation


9 Cradle

10 Pipecleaner

3 Components

No.	Classification	Description
1	Main Body	Main Body of Hand-piece
2	Tapping Rod (Applied Part)	A Tapping Rod which taps Abutment of an implant to measure implant stability and degree of Osseointegration.
3	Front Guide	Front Guide of a Tapping Rod is where the Tapping Rod contacts abutment and guides the movement of the Tapping Rod.
4	Switch	Multi-functioning Switch button that turns the device on and off and operates the tapping-motion. To turn on and off the device, press and hold the button for a while. Press the button once shortly to operate the tapping-motion.
5	LCD Display	Display function: ON / OFF display, Operation display (Check / Ready), Charging display (Charging / Charged), iST scale display (1~99), Error display (Error) *iST scale represents the degree of implant stability
6	Target Bar	Target Bar is a testing device that tests precision of measuring value. After connecting with the main body, measuring within a range of 0 ~ 30 degree is required and the number of measuring value should be in a range of ± 3 from the standard value.
7	Allen Wrench	Allen Wrench is a Wrench designed for connecting and disconnecting the Tapping Rod from the main body. While connecting and disconnecting the Tapping Rod, the Wrench enables both screwing and loosening the Rod from the main body.
8	5V/2A Adapter charger and USB C Type Cable	USB-C cable is used for charging the device by connecting the cable to the main body of the device. (If the specifications are the same, you can purchase and use a certified charging adapter)
9	Cradle	It is the place to mount the device when not in use.
10	Pipecleaner	It is used when cleaning the inside of the Front Guide.

No.	A	B	C	D
Unit (mm)	45	161.5	34.5	32
Weight (g)	100			

I. Product explanation

4 Specifications (Technical Specifications)

- ① Model Name : IMT-100
- ② Form : Internal power device & B-Type AP (Applied Part)
- ③ Electrical Rated Voltage : DC 3.7 V
- ④ Electrical Rated Current Capacitance : 640 mAh
- ⑤ Charge : USB Cable C-Type
- ⑥ Weight : 100 g
- ⑦ Size : 207mm x 34.5 mm x 32 mm
- ⑧ Operating condition : Temperature 10 °C ~ 35 °C , Humidity 35 %R.H. ~ 75 %R.H.
- ⑨ Transport and Storage condition : Temperature -20 °C ~ 40 °C , Humidity 10 %R.H. ~ 90 %R.H.
- ⑩ Operation Mode : Continuous Operation
- ⑪ Software
 - Name : IMT-SW
 - Version: IMT-1.1.X
 - Published Date : May 23rd, 2017

5 Performance

- ① Accuracy : Display value ± 3
- ② Tapping Speed : 0.5 sec. / 1 time
- ③ Tapping Strength : ≤ 1.0 N to 1.3 N
- ④ Number of Tapping : Taps 2~4 times until the measured value of 1~59. After 60 ~ 99, the device taps 6 times in order to measure the stability.

6 Symbols

- ① Safety Sign: Symbol for safety sign is displayed on exterior of the device, on package design, and on user manual. Symbol delivers important notification and caution to users. Symbols displayed below must be carefully read and followed by users prior to using the device for proper use and management.




	Warning	This Warning sign warns users that the use of damaged product may cause irreversible damage to patients.
	Caution	This Caution sign is related to possibilities of damaging product or harming patients.
	Prohibition	This Prohibition sign is related to possibilities of damaging product or harming patients.
	Mandatory Action	This Mandatory Action sign means "refer to the user manual".

② Symbol

No	Symbol	Description	Location
1		Serial No.	- Product Label - Product Package
2		Batch Code	- Product Label - Product Package
3		Refer to instruction manual/booklet	- Product Label - Product Package
4		General warning sign	- Product Package - IFU - Product label
5		Date of Manufacture	- Product Label - Product Package
6		TYPE B AP (Applied Part)	- Product Label - Product Package - IFU
7		Manufacturer	- Product Label - Product Package

I. Product explanation


6 Symbols

No	Symbol	Description	Location
8		Keep the device dry	- Product Package
9		Keep away from direct sunlight	- Product Package
10		Range of storage temperature	- Product Package
11		Range of storage humidity	- Product Package
12		Range of storage atmospheric pressure	- Product Package
13		DC Voltage	- Product Label
14		Bistable Switch	- Exterior of the device
15		Discard in the specified manner	- Product Label - Product Package - IFU
16		CE marking, Complies with european medical devices directive	- Product Label - Product Package - IFU
17		C-USA US Classification mark	- Product Label
18		Federal (USA) law restricts this device to sale by or on the order of a physician	- Product Label - Product Package
19		Indicates a medical device that has not been subjected to a sterilization process.	- Product Label - Product Package
20		Indicates the Authorized representative in the European Community.	- Product Package


I. Product explanation

6 Symbols

③ Device Label (Location: Rear of the device)


④ Target Bar Label


⑤ Package Label

⑥ Label on interior of the device

⑦ Cradle Label


II. Instructions for use

0 Frequently Used Function

- Charging with USB-C type cable
- Pressing switch in order to turn the device on and off
- Pressing Tapping Rod operation switch

1 Workflow Diagram


2 Preparation Before Use

- ① Charge at minimum 3 hours before the first use. (Battery is re-chargable)
- ② The main body of the device and the Tapping Rod must be disconnected and sterilized prior to usage. (Refer to the context below for Tapping Rod removal manual and Sterilization manual)

II. Instructions for use

3 Charging Battery (Charging adapter sold separately)

- ① Connect the USB C-Type cable, which is connected to the 5V adapter charger or 5V power terminal, to the back of the IMT-100 connection terminal.
- ② "Charging" shows up on the LCD Display while battery is charging.
- ③ "Charged" is displayed on the LCD Display when the device is fully charged.
- ④ Charging time takes approximately 3 hours depending on circumstances.

- ⑤ When the battery is fully charged, the device is capable of measuring approximately 1000 times.
- ⑥ When the battery is low, "LOW BATTERY" is displayed on the LCD Display.

- User must purchase and use appropriate charging adapters that are certified.
- Using uncertified adapters may cause a fire, explosion, or a damage to the product.
- The device is unable to use while charging.
- If the battery drains abnormally quickly or a "LOW BATTERY" sign continuously appears on the LCD Display, it indicates that the battery life is almost over. In this case, please contact our customer service for a battery exchange.
- Contact our customer service when battery does not charge for 20 minutes.
- When battery leakage occurs, stop using the device and call our customer service.

4 Operation Method

- ① Power ON
 - When the device is OFF, press the ON/OFF button once to turn on the LCD display and goes to Ready mode. (Ready mode is measuring mode).

II. Instructions for use

4 Operation Method

② Action

- When power is ON (Ready mode), shortly press the ON/OFF button to start. The device will make a short beeping sound and the tapping-motion will start operating immediately. The measured value will be displayed on the LCD Display as shown below.


③ Power OFF


- Hold the ON/OFF button for a while until the device turns off while making a short beeping sound
- If the device is not used for 2 minutes, it automatically turns off.

5 Measuring Process

- ① Turn the power on by pressing the ON/OFF switch.
- ② Carefully position the tip of the Tapping Rod to the upper edge of the healing abutment while being cautious not to push the abutment.
- ③ Maintain the contact angle within the range of $0^{\circ} \sim 30^{\circ}$ to upward while measuring.

- ④ When the Tapping Rod of the AnyCheck is located 90 degree from the side part of the top of the healing abutment and the ON/OFF switch is pressed once, the tapping motion starts after 0.5 second.
- ⑤ The measured value is more accurate when the degree between the healing abutment and the Front Guide is closer to perpendicular.
- ⑥ Check the numbers displayed on the LCD screen of the device.

- If the contact angle is off the range ($0^{\circ} \sim 30^{\circ}$ upward from the surface), an alert sound will be on and the device will not operate. In this case, tilt the device until the contact angle is back in the range; and then, try measuring again.
- Measurement error may occur if the tip of the Front Guide applies too much force to the healing abutment or is too far away from the healing abutment during measurement.

II. Instructions for use

6 Checking Measured Value (iST scale)

- ① A range of measured value is displayed as a number between 01 to 99
- ② Color of the number appears as Red for 01~59, Orange for 60~64, and Green for 65~99
- ③ Lower the measured value (iST scale), greater the implant mobility (weaker the degree of osseointegration).
- ④ The measured values differ depending on the height of the healing abutment. Calculate the correct iST scale as shown on the chart below: the standard height of the healing abutment is set as 4mm, and iST scale changes as the height of the healing abutment changes.

Standard + 3	Height 7mm Healing Abutment	Actual Value=Displayed iST scale +6
Standard + 2	Height 6mm Healing Abutment	Actual Value=Displayed iST scale +4
Standard + 1	Height 5mm Healing Abutment	Actual Value=Displayed iST scale +2
Standard	Height 4mm Healing Abutment	Actual Value=Displayed iST scale
Standard - 1	Height 3mm Healing Abutment	Actual Value=Displayed iST scale -2
Standard - 2	Height 2mm Healing Abutment	Actual Value=Displayed iST scale -4
Standard - 3	Height 1mm Healing Abutment	Actual Value=Displayed iST scale -6


- ⑤ Green light below the LCD Display is the status bar.

- The device measures iST values after each tapping up to 6 times, and if identical values are measured more than 4 times, the status bar appears green. If identical values are measured less than 4 times, the status bar appears orange. To obtain more accurate value, additional measurement is required.
- In the case that the device only taps twice (for iST value below 60), the status bar appears green if the 2 measured values are identical. On the other hand, if each measured value results unevenly, the device operates an additional tap to achieve an accurate result. Finally, if all three measured values result unevenly, the status bar displays orange.

7 Product Performance Test.

- ① Install the target bar to the Tapping Rod as shown on the picture below.

- ② After installation, measure the iST scale by pressing the switch shortly while adjusting the angle of the device into the range of upward 0° ~ 30°.
- * Error notification sound beeps when the device is trying to measure from outside of the range.

II. Instructions for use

7 Product Performance Test.

- ③ The device is functioning well if measured value on the LCD Display is (± 3) from the standard value.

* Standard value means the value written at the target bar, which is shown as same as the example below.

- ④ Please contact our customer service center if the displayed standard value shows a higher difference than (± 3) standard value more than 3 times.

8 Precautions During Use


- ① Understand the appropriate use of the product and be completely aware of the precautions prior to using the product.
- ② Do not disassemble or modify the device arbitrarily.
- ③ Do not use this product for uses other than as indicated herein.
- ④ Use at an appropriate environment within a dental clinic or hospital.
- ⑤ Use under control of a professional doctor while the patient is sitting still on the dental chair.

9 Cleaning After Use

- ① Immediately remove the Tapping Rod and the Front Guide according to the disassembling method and Immerse in ethyl alcohol (70%) and use a brush to remove foreign matter or blood.
- ② Rinse the alcohol-washed part first with water. remove any traces of blood or foreign objects remaining. Use a syringe or pipe cleaner for areas that are difficult to wash.
- ③ Final wash with flowing purified water for at least 1 minute.
- ④ Wipe it with a dry cloth or remove the water completely with warm-air circulator.(60 °C, 30 min).
- ⑤ If the surface of the body is to be cleaned, it can be wiped with a cloth soaked with water.

- Do not use disinfection sprays. It is very important to avoid the disinfection fluid penetrate into the inner parts of the product.
- Never use oil or lubricants.

10 Sterilization After cleaning

- ① In high steam sterilization(autoclave), pre-vacuum can be selected. (Only the separated "Tapping rod" and "Frongt Guide" are sterilized and reusable)
- ② The end user is to use a FDA cleared wrap for the parameter.
- ③ Sterilization condition is same below.

Cycle Type	Temperature	Pressure	Exposure Time	Dry Time
Pre-Vacuum	132 °C	2 bar	4 minutes	30 minutes
	270 °F	28.5 psi		

II. Instructions for use

10 Sterilization After cleaning


- The sterilization temperature should not exceed 135 °C (275 °F).
- Remove from autoclave immediately after finished cycle.
- Allow to cool down before use; handle carefully when hot.
- Do not wash in dishwasher.
- Do not use other cleaning fluids as they can damage the top cover of the instrument.
- Instruments are reusable for sterilization. We guarantee our products to withstand a minimum of two hundreds (200) sterilization cycles when sterilized according to the criteria listed.

11 Storage after sterilization


- ① The sterilized tapping rod must be reassembled into the IMT-100.
- ② Store in a safe place to protect against water, vibration, external shock or damages.
- ③ Store away from direct sunlight.

12 Assembling/Disassembling Method


- ① Disassemble the Front Guide from the Tapping Rod by turning the Front Guide counterclockwise as displayed below.


- ② Connect the Allen Wrench to the Tapping Rod as displayed below.


- ③ Disassemble the Tapping Rod-1 from the Tapping Rod-2 by turning the Allen Wrench counterclockwise by hand as displayed below.

- ④ For assembling the Tapping Rod, connect each parts of disassembled Tapping Rod on a reverse order.

- ⑤ Hold the Front Guide and connect it to the external metal tube by turning clockwise. Tighten the connection without any gap as displayed below.

13 Changing battery

- ① Battery must be replaced by the designated service manager.
- ② The specifications of the battery to be replaced are as follows.
 - Manufacturer: Hangzhou Future Power Technology Co., Ltd
 - Type: Lithium polymer
 - Model Name: FT622031P
 - Rating: 3.7 Vdc
 - Volume: 640 mAh
- ③ Follow the process mentioned earlier at number 12 (Assembling/ Disassembling Method) to remove the Tapping Rod.
- ④ Disassemble the exterior by sliding the upper outer part and the lower outer part of the main body.
- ⑤ Disassemble the connector from the battery that needs to be replaced.
- ⑥ Check the specifications of the new battery; and then, connect the new battery to the connector. Close the battery cover by sliding the cover other way around.
- ⑦ Follow the process mentioned earlier at number 12 to assemble the Tapping Rod.
- ⑧ The device performance must be tested after battery replacement.

14 Storage and Maintenance After Use (Instructions to stop operation)

- ① Turn off the device by press and holding the ON/OFF switch. Then, disassemble the Applied Part(AP).
- ② Clean the Hand Piece part of the device with alcohol cloth and charge the device by connecting it to the USB C Type cable.
- ③ Clean and Store the AP.
- ④ Clean the Tapping Rod and Front Guide in contact with the patient after use, and sterilize them before use. For cleaning and sterilization of the Tapping Rod and Front Guide, assembly and disassembly methods shall be carried out in accordance with the instructions for use. Cleaning process for the device is same as the process for regular medical instruments. and sterilization is performed according to the sterilization conditions presented in the manual

15 Disposal of Products and Components

- Device and USB cable shall be disposed by following the Disposal Guideline or returned to the manufacturer.
- Charging adapter must be disposed by following the disposal process suggested by the supplier.

III. Warnings and Safety Precautions

1 Precaution

- ① The device stops operating and makes error notification sound when measuring angle is off the range of 0 ~ 30 degree upward based on the ground level.
- ② The tip of the Tapping Rod is designed to contact the healing abutment for the convenience of users. Applying too much pressure onto the healing abutment or holding the device too far away from the healing abutment during measurement may result in an error.
- ③ Device must be charged when the ON/OFF switch is not functioning. (* Device does not work while charging)
- ④ Battery has to be replaced to new battery when the device turns off abnormally or "LOW BATTERY" sign continuously appear on the LCD screen.
- ⑤ User must contact the customer service center when the battery needs to be replaced since the battery is built-in to the device.
- ⑥ Must be used within dental hospital.
- ⑦ Use under control of a professional doctor while the patient is sitting still on the dental chair.

2 Prohibition

- ① Non-professionals are prohibited to use the device.
- ② It should not be used for any purpose other than diagnosis for implant treatment
- ③ Prohibited to replace any parts of the device for recreation or to connect to anything other than the parts supplied by the manufacturer.
- ④ Prohibited to use near flammable or dangerous materials.
- ⑤ Prohibited to use the product in an environment where electromagnetic waves are exposed as this may cause malfunction.
- ⑥ Prohibited to use near water as the device is not waterproof.

3 Warnings

- ① Be sure to fully understand the operating instructions and precautions for the product before use.
- ② An authorized charging adapter should be used for the sake of safety.
- ③ The product should be inspected by the manufacturer if the product is impacted.
- ④ The device must be carried within the container offered with the package.
- ⑤ The product is not disinfected or sterilized at the time of initial delivery; thus, the device must be used after sterilization.
- ⑥ In patients with high cortical bone thickness and bone density, several factors related to the stability of the implants should be considered.
- ⑦ Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- ⑧ Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

IV. Statements and tables for EMC

Table 1 - ELECTROMAGNETIC EMISSIONS – for IMT-100

Guidance and manufacturer's declaration – electromagnetic emissions

IMT-100 is intended for use in the electromagnetic environment specified below. The customer or the user of the IMT-100 should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The IMT-100 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The IMT-100 is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded: Warning: This equipment/system is intended for use by professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the IMT-100 or shielding the location.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 2 – ELECTROMAGNETIC IMMUNITY – for IMT-100

Guidance and manufacturer's declaration – electromagnetic immunity

The IMT-100 is intended for use in the electromagnetic environment specified below. The customer or the user of the IMT-100 should assure that it is used in such an environment.

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) EN 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst EN 61000-4-4	N/A	N/A	N/A
Surge EN 61000-4-5	N/A	N/A	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines EN 61000-4-11	N/A	N/A	N/A
Power frequency (50/60 Hz) magnetic field EN 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE UT is the a.c. mains voltage prior to application of the test level.

IV. Statements and tables for EMC

Table 3 - ELECTROMAGNETIC – for IMT-100 that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration – electromagnetic immunity

The IMT-100 is intended for use in the electromagnetic environment specified below. The customer or the user of the IMT-100 should assure that it is used in such an environment.

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF EN 61000-4-6	N/A	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the [ME EQUIPMENT or ME SYSTEM], including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $\sqrt{P}$ $d = 1.17$ $\sqrt{P}$ $d = 1.17 \text{ 80 MHz to 800 MHz}$ $\sqrt{P}$ $d = 2.33 \text{ 80 MHz to 800 MHz}$ where P is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer and d is the recommended separation distance in metres(m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF EN 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the IMT-100 is used exceeds the applicable RF compliance level above, the IMT-100 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the IMT-100.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

IV. Statements and tables for EMC

Table 4 – Recommended separation distances between portable and mobile RF communications equipment and the ME EQUIPMENT or ME SYSTEM – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the IMT-100

The IMT-100 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the IMT-100 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the IMT-100 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $\sqrt{P}$ d = 1.17	80 MHz to 800 MHz $\sqrt{P}$ d = 1.17	800 MHz to 2,5 GHz $\sqrt{P}$ d = 2.33
0.01	0.117	0.117	0.233
0.1	0.370	0.370	0.736
1	1.17	1.17	2.33
10	3.70	3.70	7.36
100	11.7	11.7	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

V. Statements for FCC

1. Certification and Safety Approvals

1.1 FCC Compliance Statement

This device complies with part 15 of the FCC rules.

Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment on and off, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antennae
- Increase the separation between the equipment and the receiver
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

1.2 FCC RF Exposure Statement

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. End users must follow the specific operating instructions for satisfying RF exposure compliance. The antenna used for this transmitter must not transmit simultaneously with any other antenna or transmitter, except in accordance with FCC multi-transmitter product procedures.

1.3 FCC Caution

Any changes or modifications to the equipment not expressly approved by the party responsible for compliance could void user's authority to operate the equipment.

1.4 Industry Canada Statement

This device complies with Industry Canada license-exempt RSS standard(s).

Operation is subject to the following two conditions:

- (1) This device may not cause interference.
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

1.5 This device can be operated in at least one MemberState without infringing applicable requirements on the use of radio spectrum.

- RF Range: 2402 ~ 2480MHz

- Max Output Power(E.I.R.P.): "ant gain + RF output power" dBm

VI. Product Warranty Policy


- This product is manufactured under strict inspection and thorough product management by DMS Co., Ltd.
- Compensation standards for product repairs and exchanges are followed by the "Consumer Injury Compensation Rules"
- DMS Co., Ltd. warrants that reasonable care has been used in the design and manufacture of this product. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness.
- Handling, storage and cleaning of this product as well as factors relating to the patient, diagnosis and other matters beyond DMS Co., Ltd.'s control directly affect the product and the results obtained from its use.
- DMS Co., Ltd.'s obligation under this warranty is limited to the repair or replacement of this product and DMS Co., Ltd. shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this product.
- DMS Co., Ltd. neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this product. DMS Co., Ltd. assumes no liability with respect to products reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for intended use, with respect to such product.
- Contact the address below for more information about the product service
Manufacturer: DMS Co., Ltd.
Address: 4dong-#107, 201, 202, 42-10, Taejunggongdan-gil, Wonju-si, Gangwon-do, Republic of Korea
Tel: +82-33-742-3210 Fax: +82-33-742-3219

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- Este producto se ha fabricado bajo la estricta revisión y el exhaustivo control de productos de DMS Co., Ltd.
Fabricante: DMS Co., Ltd.